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1. National Institute of Environmental Health Sciences. NIEHS Perfluorinated Chemicals (PFCs) fact sheet, 2012. https://www.atsdr.cdc.gov/sites/pease/documents/perfluorinated_chemicals_508.pdf (accessed Dec 1, 2019).
2. United States Environmental Protection Agency. Per- and Polyfluoroalkyl Substances (PFAS). <https://www.epa.gov/pfas> (accessed Dec 1, 2019).
3. Agency for Toxic Substances and Disease Registry (ATSDR). An Overview of Perfluoroalkyl and Polyfluoroalkyl Substances and Interim Guidance for Clinicians Responding to Patient Exposure Concerns, 2019. https://www.atsdr.cdc.gov/pfas/docs/ATSDR_PFAS_ClinicalGuidance_12202019.pdf (accessed Apr 6, 2020).
4. DIN e.V. DIN 38407-42:2011-03 - Jointly determinable substances (group F) - Part 42: Determination of selected polyfluorinated compounds (PFC) in water - Method using high performance liquid chromatography and mass spectrometric detection (HPLC/MS-MS) after SPE; Beuth Verlag: Berlin, 2011.
5. U.S Environmental Protection Agency (EPA). Method 537_FINAL_rev1.1 Determination of Selected Perfluorinated Alkyl Acids in Drinking Water by Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS) (Sept,2009), 2009. <https://nepis.epa.gov/Exe/ZyPDF.cgi?Dockey=P100EQ6W.txt> (accessed December 1, 2019).
6. ASTM International. ASTM D7979-17, Standard Test Method for Determination of Per- and Polyfluoroalkyl; Conshohocken, PA, 2017.
7. Organisation for Economic Co-operation and Development. Updated List complete document 25-06-07 - revised August 20, 2007. Lists of PFOS, PFAS, PFOA, PFCA, Related Compounds and Chemicals that may degrade to PFCA. [http://www.oecd.org/officialdocuments/publicdisplaydocumentpdf/?doclanguage=en&cote=en/jm/mono\(2006\)15](http://www.oecd.org/officialdocuments/publicdisplaydocumentpdf/?doclanguage=en&cote=en/jm/mono(2006)15) (accessed Dec 1, 2019).
8. Organisation for Economic Co-operation and Development. Toward a new comprehensive global database of per- and polyfluoroalkyl substances (PFASs): Summary report on updating the OECD 2007 Lists of Per- and Polyfluoroalkyl substances (PFASs), 2018. [https://www.oecd.org/officialdocuments/publicdisplaydocumentpdf/?cote=ENV-JM-MONO\(2018\)7&doclanguage=en](https://www.oecd.org/officialdocuments/publicdisplaydocumentpdf/?cote=ENV-JM-MONO(2018)7&doclanguage=en) (accessed Jan 6, 2020).
9. Houtz, E. F.; Sutton, R.; Park, J.-S.; Sedlak, M. Poly- and perfluoroalkyl substances in wastewater: Significance of unknown precursors, manufacturing shifts, and likely AFFF impacts. *Water Res.*, **2016**, *95*, 142–149.
10. D'Agostino, L. A.; Mabury, S. A. Certain Perfluoroalkyl and Polyfluoroalkyl Substances Associated with Aqueous Film Forming Foam Are Widespread in

- Canadian Surface Waters. *Environ. Sci. Technol.*, **2017**, *51*, 13603–13613.
11. Abwasserverordnung in der Fassung der Bekanntmachung vom 17. Juni 2004 (BGBl. I S. 1108, 2625), die zuletzt durch Artikel 1 der Verordnung vom 22. August 2018 (BGBl. I S. 1327) geändert worden ist.
 12. Oberflächengewässerverordnung vom 20. Juni 2016 (BGBl. I S. 1373).
 13. Hu, J.; Rohrer, J. Determination of Adsorbable Organic Halogen in Wastewater Using a Combustion Ion Chromatography System, 2017. Application Note 72333. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-72333-IC-Adsorbable-Organic-Halogen-Wastewater-AN72333-EN.pdf> (accessed Dec 1, 2019).
 14. Yang, H.; Khor, D.; Rohrer, J. Determination of fluoride in tea using a combustion ion chromatography system, 2017. Application Note 72268. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-72268-IC-Fluoride-Tea-AN72268-EN.pdf> (accessed Dec 12, 2019).
 15. Aggrawal, M.; Rohrer, J. Determination of Chlorine, Bromine, and Sulfur in Polyethylene Materials Using Combustion Ion Chromatography, 2017. Application Note 72349. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-72349-IC-Chlorine-Bromine-Sulfur-Polyethylene-AN72349-EN.pdf> (accessed Dec 12, 2019).
 16. Agustin, A.; Fisher, C.; Manahan, M. Determination of Halogens in Coal Using Combustion Ion Chromatography, 2017. Application Note 1145. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-1145-IC-Halogens-Coal-AN1145-EN.pdf> (accessed Dec 12, 2019).
 17. Aggrawal, M.; Rohrer, J. Determination of Chlorine, Bromine, and Sulfur in Polyethylene Materials Using Combustion IC with a Carbonate/Bicarbonate Eluent, 2018. Application Update 72588. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/au-72588-ic-combustion-cl-br-s-polyethylene-au72588-en.pdf> (accessed Dec 12, 2019).
 18. International Organization for Standardization. *ISO 10304-1:2007-08 - Water quality - Determination of dissolved anions by liquid chromatography of ions - Part 1: Determination of bromide, chloride, fluoride, nitrate, nitrite, phosphate and sulfate*; Beuth Verlag: Berlin, 2007.
 19. U.S. Environmental Protection Agency. Method 1650: Adsorbable Organic Halides by Adsorption and Coulometric Titration; Revision C, 1997. https://www.epa.gov/sites/production/files/2015-10/documents/method_1650c_1997.pdf (accessed December 1, 2019).
 20. US Environmental Protection Agency. Method 9020B: Total Organic Halides (TOX), part of Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, 1994. <https://www.epa.gov/sites/production/files/2015-12/documents/9020b.pdf> (accessed Dec 1, 2019).
 21. DIN e.V. *Working Draft DIN 38409-59, to be published. Details to be obtained from Franz Schmitz, mail to: franzschmitz.wi@t-online.de.*
 22. von Abercron, E.; Falk, S.; Stahl, T.; Georgii, S.; Hamscher, G.; Brunn, H.; Schmitz, F. Determination of adsorbable organically bound fluorine (AOF) and adsorbable organically bound halogens as sum parameters in aqueous environmental samples using combustion ion chromatography (CIC). *Sci. Total Environ.*, **2019**, *673*, 381–391.
 23. Mitsubishi Chemical Analytech Co., Ltd. AQF-2100H Automatic Quick Furnace Combustion Ion Chromatography, 2014. AQF 2100H Brochure - COSA Xentaur. https://cosaxentaur.com/resources/files/1038/AQF_2100H_COSA_Xentaur.pdf (accessed Jan 13, 2020).
 24. Khor, D.; Rohrer, J.; Agustin, A.; Manahan, M. Combustion Ion Chromatography with a Dionex Integrion HPLC System, 2019. Technical Note 72211. <https://assets.thermofisher.com/TFS-Assets/CMD/Technical-Notes/tn-72211-integrion-hpic-combustion-ic-tn72211-en.pdf> (accessed Dec 12, 2019).
 25. International Organization for Standardization. *ISO 9562:2004-09-29 - Water quality - Determination of adsorbable organically bound halogens (AOX)*; Beuth-Verlag: Berlin, 2004.
 26. DIN e.V. *DIN 38402-51:2017-05 - German standard methods for the examination of water, waste water and sludge - General information (group A) - Part 51: Calibration of analytical methods - Linear calibration (A 51)*; Beuth Verlag: Berlin, 2017.
 27. DIN e.V. *DIN 38402-60:2013-12 - German standard methods for the examination of water, waste water and sludge - General information (group A) - Part 60: Analytical quality assurance for chemical and physicochemical water analysis (A 60)*; Beuth Verlag: Berlin, 2013.
 28. Wagner, A.; Raue, B.; Brauch, H.-J.; Worch, E.; Lange, F. T. Determination of adsorbable organic fluorine from aqueous environmental samples by adsorption to polystyrene-divinylbenzene based activated carbon and combustion ion chromatography. *J. Chromatogr. A*, **2013**, *1295*, 82–89.
 29. Willach, S.; Brauch, H.-J.; Lange, F. T. Contribution of selected perfluoroalkyl and polyfluoroalkyl substances to the adsorbable organically bound fluorine in German rivers and in a to the adsorbable organically bound fluorine in German rivers and in a highly contaminated groundwater. *Chemosphere*, **2016**, *145*, 342–350.
 30. Furl, C. V.; Meredith, C. A.; Strynar, M. J.; Nakayama, S. F. Relative importance of wastewater treatment plants and non-point sources of perfluorinated compounds to Washington State rivers. *Sci. Total Environ.* **2011**, *409*(15), 2902–2907.
 31. Murakami, M.; Shinohara, H.; Takada, H. Evaluation of wastewater and street runoff as sources of perfluorinated surfactants (PFSS). *Chemosphere*, **2009**, *74*(4), 487–493.
 32. Zushi, Y.; Takeda, T.; Masunaga, S. Existence of nonpoint source of perfluorinated compounds and their loads in the Tsurumi River basin, Japan. *Chemosphere*, **2008**, *71*(8), 1566–1573.
 33. Zushi, Y.; Masunaga, S. Identifying the nonpoint source of perfluorinated compounds using a geographic information system based approach. *Environ. Toxicol. Chem.*, **2009**, *28*(4), 691–700.
 34. Houtz, E. F.; Sedlak, D. L. Oxidative Conversion as a Means of Detecting Precursors to Perfluoroalkyl Acids in Urban Runoff. *Environ. Sci. Technol.*, **2012**, *46*, 9342–9349.

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© 2021 Thermo Fisher Scientific Inc. All rights reserved.
All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.