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1. Applications of Ion Chromatography for Pharmaceutical and Biological Products (ed.: L. Bhattacharyya and J. S. Rohrer), John Wiley and Sons Inc. (Hoboken, New Jersey) (2012).
2. Thermo Scientific Application Note 148: Determination of bethanechol by ion chromatography. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-148-IC-Bethanechol-LPN1510-EN.pdf> (accessed Mar. 25, 2022)
3. Li, Y.; Nesterenko, P.; Stanley, R.; Paull, B.; Macka, M. Comparison of cation-exchange capillary columns used for ion chromatographic separation of biogenic amines. *Journal of Chromatography A* **2018**, *1571*, 193–200.
4. Thermo Scientific Dionex IonPac IC column selection guide. <https://assets.thermofisher.com/TFS-Assets/CMD/brochures/st-70588-ic-column-selection-guide-st70588-en.pdf> (accessed Mar. 25, 2022)
5. Thermo Scientific Application Note 1073: Determination of ammonia in sodium bicarbonate. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-1073-IC-Ammonia-Sodium-Bicarbonate-AN70815-EN.pdf> (accessed Mar. 25, 2022)
6. Thermo Scientific Application Note 73482: Determination of ammonia impurity in potassium bitartrate using ion chromatography. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/an-73482-ic-ammonia-impurity-potassium-bitartrate-an73482-en.pdf> (accessed Mar. 25, 2022)
7. Thermo Scientific Application Note 194: Determination of carbachol in ophthalmic solutions using a Reagent-Free ion chromatography system. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-194-IC-Carbachol-Ophthalmic-Solutions-LPN1967-EN.pdf> (accessed Mar. 25, 2022)
8. Thermo Scientific Application Note 199: Determination of N-methylpyrrolidine in cefepime Using a Reagent-Free ion chromatography system. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-199-IC-N-Methylpyrrolidine-Cefepime-LPN2005-EN.pdf> (accessed Mar. 25, 2022)
9. Thermo Scientific Application Note 249: Determination of methacholine chloride and potential impurities using a Reagent-Free ion chromatography system. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-249-IC-Methacholine-Chloride-LPN2517-EN.pdf> (accessed Mar. 25, 2022)
10. Thermo Scientific Application Note 298: Determination of dimethylamine in metformin HCl drug product using IC with suppressed conductivity detection. <https://appslab.thermofisher.com/App/1474/an298-determination-dimethylamine-metformin-hcl-drug-product-using-ic-with-suppressed-conductivity-detection> (accessed Mar. 25, 2022)
11. Thermo Scientific Application Note 1057: Determination of methylamine in drug products. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-1057-IC-Monomethylamine-Drugs-AN70523-EN.pdf> (accessed Mar. 25, 2022)

12. Thermo Scientific Application Note 1062: Determination of morpholine in linezolid by ion chromatography. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-1062-IC-Morpholine-Linezolid-AN70702-EN.pdf> (accessed Mar. 25, 2022)
13. Thermo Scientific Application Note 72649: Validation of an improved ion chromatography method for the limit of choline test in the USP Succinylcholine Chloride monograph. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/an-72649-ic-limit-choline-test-an72649-en.pdf> (accessed Mar. 25, 2022)
14. Concise International Chemical Assessment Document, N-NITROSODIMETHYLAMINE. <https://www.who.int/ipcs/publications/cicad/en/cicad38.pdf?ua=1> (accessed Mar. 25, 2022)
15. U.S. Environmental Protection Agency, Integrated Risk Information System (IRIS). https://cfpub.epa.gov/ncea/iris/iris_documents/documents/subst/0045_summary.pdf (accessed Mar. 25, 2022)
16. FDA announces voluntary recall of several medicines containing valsartan following detection of an impurity. <https://www.fda.gov/news-events/press-announcements/fda-announces-voluntary-recall-several-medicines-containing-valsartan-following-detection-impurity> (accessed Mar. 25, 2022)
17. FDA Updates and Press Announcements on NDMA in Zantac (ranitidine). <https://www.fda.gov/drugs/drug-safety-and-availability/fda-updates-and-press-announcements-ndma-zantac-ranitidine> (accessed Mar. 25, 2022)
18. Mylan initiates voluntary nationwide recall of three lots of nizatidine capsules, USP, due to the detection of trace amounts of NDMA (N-nitrosodimethylamine) impurity found in the active pharmaceutical ingredient manufactured by Solara Active Pharma Sciences Limited. <https://www.fda.gov/safety/recalls-market-withdrawals-safety-alerts/mylan-initiates-voluntary-nationwide-recall-three-lots-nizatidine-capsules-usp-due-detection-trace> (accessed Mar. 25, 2022)
19. FDA Updates and Press Announcements on NDMA in Metformin. <https://www.fda.gov/drugs/drug-safety-and-availability/fda-updates-and-press-announcements-ndma-metformin> (accessed Mar. 25, 2022)
20. GUIDANCE DOCUMENT Control of Nitrosamine Impurities in Human Drugs Guidance for Industry. <https://www.fda.gov/regulatory-information/search-fda-guidance-documents/control-nitrosamine-impurities-human-drugs> (accessed Mar. 25, 2022)
21. USP General Chapter Prospectus: <1469> Nitrosamine Impurities. <https://www.uspnf.com/notices/nitrosamine-impurities-gc-prospectus-20200424> (accessed Mar. 25, 2022)
22. Mirvish, S. S. Formation of N-nitroso compounds: Chemistry, kinetics and in vivo occurrence. *Toxicol. Appl. Pharmacol.* **1975**, 31, 325–351.
23. Ashworth, I. W.; Dirat, O.; Teasdale, A.; Whiting, M. Potential for the formation of N-nitrosamines during the manufacture of active pharmaceutical ingredients: an assessment of the risk posed by trace nitrite in water. *Org. Process Res. Dev.* **2020**, DOI: 10.1021/acs.oprd.0c00224
24. Thermo Scientific Application Note 73987: 医薬品における亜硝酸イオンの測定 <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/an-73987-ic-nitrite-pharmaceuticals-an73987-ja.pdf>
25. Thermo Scientific Dionex ICS-5000+ Ion Chromatography System Operator's Manual December 2014. <https://assets.thermofisher.com/TFS-Assets/CMD/manuals/man-065446-ics-5000-plus-man065446-en.pdf> (accessed Mar. 25, 2022)
26. Thermo Scientific Dionex IonPac CS19 column product manual. <https://assets.thermofisher.com/TFS-Assets/CMD/manuals/man-065440-ionpac-cs19-columns-man065440-en.pdf> (accessed Mar. 25, 2022)
27. ICH Guideline Q2B, Validation of Analytical Procedures, Methodology (CPMP/ICH/281/95), Geneva, Switzerland, November 1996.

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TEL : 0120-753-670 FAX : 0120-753-671
 Analyze.jp@thermofisher.com

facebook.com/ThermoFisherJapan @ThermoFisherJP

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© 2022 Thermo Fisher Scientific Inc. All rights reserved.
All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.
Sigma-Aldrich is a trademark of Merck KGaA.
Nytol and QuickCap are trademarks of Medtech Products Inc.
Benadryl is a trademark of Johnson & Johnson Corporation.